



Quality of Life Improvement



Prevent Stress Urinary Incontinence



Post-partum Recovery



Rehabilitation of Pelvic Muscles



Pubococcygeous Enhancement



Women's Private Health Improvement

UMSELLA[™]

A breakthrough therapy for incontinence
A revolution technology for women's private health

● Tel: 0086 - 371 - 55677868 0086 - 371 - 55677886 0086 - 371 - 55011989

Headquarters:

Building 2, Jiatu Property, Hongsong Road, Lianhua Street, High-tech Zone, Zhengzhou, China

For more information: www.pzlaser.com E-mail: info@pzlaser.com

UMSELLA™



Smart On/Off System



Entirely Non-Invasive



Remarkable Effect



No Pain, No Downtime




Clinically Proven



WHAT IS UMSELLA?

UMSELLA uses high-intensity focused electromagnetic technology to stimulate all the deep muscles of the pelvic floor. During the brief treatment, muscles contract in response to the electromagnetic stimulation, producing the same result as doing 12,000 kegel exercises.



BEFORE

Pelvic floor muscles insufficiently support pelvic organs and affect bladder control.

TREATMENT

UMSELLA™ effectively stimulates pelvic floor muscles

AFTER

Stimulation leads to regained control over pelvic floor muscles and bladder.

WHY UMSELLA IS UNIQUE?

Smart On/Off System

With ergonomic benefits, the slide-out technology will reduce risk of overreaching - there when you need it, concealed when you don't.

Entirely Non-Invasive

While remaining fully clothed, women sit on the UMSELLA chair for about 30 minutes. It sounds simple, and it is: no surgery, no medication and completely non-invasive.

Remarkable Effect

You may observe improvement after a single session. The results will typically continue to improve over the next few weeks.

No Pain, No Downtime

You will experience tingling and pelvic floor muscle contractions during the procedure. You may resume daily activities immediately after the treatment.


Clinically Proven

Scientific research has shown that 95% of treated patients reported significant improvement in their quality of life.




RESULTS

The graph below shows the results of therapy using functional magnetic stimulation after 16 treatment sessions. The results were collected in three independent centers with 160 female patients.



MEASUREMENTS OF URODYNAMICS were performed on a female patient before and after 16 treatment sessions.



Legend:
 — Measurement of the pressure within the bladder, before using the UMSELLA
 — Measurement of the pressure within the bladder, after using the UMSELLA

WHO IS THE RIGHT CANDIDATE?

UMSELLA is a great option for women of any age who desire solution for urinary incontinence and improvement in their quality of life.

Body aging, childbirth, and menopause lead to incontinence

Stress Urinary Incontinence (SUI)

Overactive Bladder

Prostatic hyperplasia

Mild to moderate bladder leakage

Postpartum recovery

Pelvic floor muscle improvement

Quality of life enhancement

TECHNICAL ADTA:

UMSELLA™	
Voltage	200-240V/100V-120V 50/60HZ
Pulse Width	200us
Power	1800W
Waveform	biphasic wave
Frenquency	1-150HZ

